

EXREsolve™ Prestained Protein Markers 8-250 kDa

Catalog Number: EXBR018



info@exreprotein.com

For Research Use Only. Not Intended for Diagnostic or Therapeutic Use.

Product Description

The EXREsolve™ Prestained Protein Markers product is composed of 11 purified pre-stained proteins with a molecular weight range of 8 kDa to 250 kDa. The individual protein markers are 8 kDa, 17 kDa, 25 kDa, 33 kDa, 43 kDa, 55 kDa, 72 kDa, 100 kDa, 130 kDa, 180 kDa, and 250 kDa. The 8 kDa marker is green, the 72 kDa marker is orange, and the other 9 markers are blue. This product is suitable as a protein molecular weight standard for SDS-PAGE or Western blot and is compatible with PVDF, nylon, and nitrocellulose membranes. The prestained markers enable direct observation of the protein electrophoresis status and estimation of protein transfer during Western Blot.

The product is a ready-to-use liquid that does not require heating, dilution or additional reducing agents before loading sample. The concentration of each protein band in the product is approximately 0.1-0.4 mg/mL.

Available sizes: 1 Pack (250 uL) or 5 Pack (5 x 250 uL)

For the separation of 10-250 kDa proteins, we recommend EXREsolve™ FlexKD PAGE Gel (catalog # EXBR001).

Limitations

- For Research Use Only. Not Intended for Diagnostic or Therapeutic Use.

Precautions

- Always wear appropriate protective clothing and follow safe laboratory procedures.

Storage

- Shipped on blue ice.
- Store at -20°C for up to 12 months. **Avoid repeated freeze-thaw cycles.**
- **Storage Buffer:** 62.5 mM Tris-H₃PO₄ (pH 7.5 at 25°C), 5 mM EDTA, 1% (w/v) SDS, 1 mM DTT, 0.02% (v/v) ProClin™300, 33% (v/v) glycerol.

Protocol

1. Thaw at room temperature for a few minutes prior to use. **Do not boil.**
2. Load 5 uL per well on an SDS-polyacrylamide gel.
 - Volumes may need to be increased for low-percentage gels (<8%), large gels, or gel thickness ≥ 1.5 mm.

- Note:**
- Transfer times of high molecular weight proteins (>100 kDa) during Western blotting may need to be increased.
 - The low molecular weight protein markers may migrate with the dye front in low-percentage gels (<8%).
 - The mobility of the prestained protein markers may be affected by various buffers and gel concentrations. To achieve suitable approximation of molecular weight determination, it is recommended to calibrate against unstained standards using the same conditions.

Figure 1. EXREsolve™ Prestained Protein Markers 8-250 kDa Typical Data




Figure 2. Example of EXREsolve™ Prestained Protein Markers 8-250 kDa Mobility in Various Conditions

Gel type	Tris-Glycine						Bis-Tris						Tris-Acetate		Hepes-Tris		
	8%	10%	12.5%	15%	B4-20%	W4-20%	G4-12%	G8-16%	G4-20%	G4-12%	G8-16%	G4-20%	G10%	6%	T3-8%	W4-20%	
Running buffer	Tris-Glycine						MES			MOPS			Tris-Acetate		Hepes		
Apparent Molecular Weights, kDa																	
% length of gel	10	250	250	250	250	250	250	250	250	250	250	250	250	250	250	250	250
	20	180	130	130	100	250	180	180	180	180	180	180	180	180	180	180	180
	30	130	72	75	55	180	130	130	130	130	130	130	130	130	130	130	130
	40	100	55	43	33	130	100	100	100	100	100	100	100	100	100	100	100
	50	70	43	33	25	100	70	70	70	70	70	70	70	70	70	70	70
	60	33	33	33	25	70	55	55	55	55	55	55	55	55	55	55	55
	70	55	25	25	17	55	43	43	43	43	43	43	43	43	43	43	43
	80	43	17	17	17	43	33	33	33	33	33	33	33	33	33	33	33
	90	17	8	8	8	33	25	25	25	25	25	25	25	25	25	25	25
	100	33	8	8	8	25	17	17	17	17	17	17	17	17	17	17	17

ProClin is a trademark of The Dow Chemical Company or an affiliated company of Dow.